



# The Protector Line

## SAFE- - SHIELD



[www.trianglerubberco.com](http://www.trianglerubberco.com)

### Rubber Visor

(Latex free polyisoprene / EPDM rubber)


Adjustable rubber straps

Rubber grommets to attach shield

### Clear Shield

(4-6 mil optical film, low distortion, optically clear)

This item is replaceable up to 5 times



**Triangle Rubber and Plastics has been serving the Automotive, Appliance, and Health and Safety industries for over 65 years.**

**We are committed to helping our community with this new line of safety products.**

This is a face shield to be used for medical purposes. This product is not integrated with any other PPE such as a facemask, but rather is for use as a standalone face shield. This product is intended for a single user. The visor can be sanitized per CDC or EPA sanitization recommendations for PPE face shields. The clear shields are intended for a single use. It is recommended that they be removed, the visor is sanitized and a new shield be assembled to the visor. Not recommended for use in the presence of high intensity heat source or flammable gas. This product has not been FDA Approved. The product has been authorized by FDA under an EUA for use by healthcare providers as personal protective equipment. This product is only authorized for the duration of the declaration that circumstances justifying the authorization of emergency use under section 564(b)(1) of the Act, 21 USC 360bbb-3(b)(1).